

产品名称: PI3 Kinase p55 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83778



产品概述 (Summary)

产品名称 (Production Name)	PI3 Kinase p55 gamma Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PI3 Kinase p55 gamma
别名 (Alternative Names)	GRB1, P85A, PI3-kinase p85-alpha subunit, PI3K, PI3K p85-alpha; P55G, PI3-kinase p85-gamma subunit, PI3K p85-gamma,;PI3 kinase p55 gamma
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92569.A synthesized peptide derived from human PI3 kinase p55 gamma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	54 kDa

产品名称: PI3 Kinase p55 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83778




研究背景 (Background)


Binds to activated (phosphorylated) protein-tyrosine kinases through its SH2 domain and regulates their kinase activity. During insulin stimulation, it also binds to IRS-1.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




产品名称: PI3 Kinase p55 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83778




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunofluorescent analysis of C6 cells, using PI 3 Kinase p55 gamma Antibody .




Immunohistochemical analysis of paraffin-embedded Human pituitary tumor, using the Antibody at 1:200 dilution.




Immunohistochemical analysis of paraffin-embedded Human lung adenocarcinoma, using the Antibody at 1:200 dilution.


产品名称: PI3 Kinase p55 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83778



Immunohistochemical analysis of paraffin-embedded Human pituitary tumor, using the Antibody at 1:200 dilution.



Immunohistochemical analysis of paraffin-embedded Human liver cancer, using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded human astrocytoma, using PI 3 Kinase p55 gamma Antibody.

注意事项 (Note)

For research use only .