

产品名称: DGCR8 Rabbit Monoclonal Antibody
产品货号: AMRe83786



产品概述 (Summary)

产品名称 (Production Name)	DGCR8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DGCR8
别名 (Alternative Names)	DGCR8; DGCRK6; Gy1 ; C22orf12; D16Wis2; pasha;; DGCR8
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8WYQ5.A synthesized peptide derived from human DGCR8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 86 kDa ; Observed MW: 100 kDa

研究背景 (Background)

Component of the microprocessor complex that acts as a RNA- and heme-binding protein that is involved in the initial step


产品名称: DGCR8 Rabbit Monoclonal Antibody
产品货号: AMRe83786




of microRNA (miRNA) biogenesis. Component of the microprocessor complex that is required to process primary miRNA transcripts (pri-miRNAs) to release precursor miRNA (pre-miRNA) in the nucleus.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.


产品名称: DGCR8 Rabbit Monoclonal Antibody
产品货号: AMRe83786



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .