

产品名称: HIF-1beta Rabbit Monoclonal Antibody
产品货号: AMRe83787



产品概述 (Summary)

产品名称 (Production Name)	HIF-1beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HIF-1beta
别名 (Alternative Names)	ARNT;HIF-1beta;HIF1B;HIF1BETA;TANGO;bHLHe2;;HIF 1 beta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P27540.A synthesized peptide derived from human HIF 1 beta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	87 kDa

研究背景 (Background)

Required for activity of the AHR. Upon ligand binding, AHR translocates into the nucleus, where it heterodimerizes with


产品名称: HIF-1beta Rabbit Monoclonal Antibody
产品货号: AMRe83787




ARNT and induces transcription by binding to xenobiotic response elements (XRE). Not required for the ligand-binding subunit to translocate from the cytosol to the nucleus after ligand binding.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: HIF-1beta Rabbit Monoclonal Antibody
产品货号: AMRe83787



Immunohistochemical analysis of paraffin-embedded human colon, using HIF-1 beta Antibody .

注意事项 (Note)

For research use only .