

产品名称: Phospho-JNK1/2/3(T183+T183+T221) Rabbit
Monoclonal Antibody
产品货号: AMRe83788



产品概述 (Summary)

产品名称 (Production Name)	Phospho-JNK1/2/3(T183+T183+T221) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-JNK1/2/3(T183+T183+T221)
别名 (Alternative Names)	JNK 46; JNK 55; MAPK10; MAPK9; MAPK8; SAPK1b ; SAPK1; SAPK; PRKM10; PRKM9; PRKM8;;p-JNK1/2/3 (T183/T183/T221)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P45983/P45984/P53779.A synthesized peptide derived from human JNK1 around the phosphorylation site of T183

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
----------------------	--

产品名称: Phospho-JNK1/2/3(T183+T183+T221) Rabbit
Monoclonal Antibody
产品货号: AMRe83788


蛋白分子量 (Molecular Weight) Calculated MW: 48,53 kDa ; Observed MW: 46,54 kDa

研究背景 (Background)


Serine/threonine-protein kinase involved in various processes such as cell proliferation, differentiation, migration, transformation and programmed cell death. Extracellular stimuli such as proinflammatory cytokines or physical stress stimulate the stress-activated protein kinase/c-Jun N-terminal kinase (SAP/JNK) signaling pathway.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.


产品名称: Phospho-JNK1/2/3(T183+T183+T221) Rabbit
Monoclonal Antibody
产品货号: AMRe83788



Western blot analysis of JNK1/2/3 phosphorylation expression in NIH/3T3 cell lysate treated with Anisomycin.




Immunofluorescent analysis using the Antibody at 1:50 dilution.




Immunofluorescent analysis of NIH/3T3 cells treated with Anisomycin, using Phospho-JNK1/2/3 (T183+T183+T221)
Antibody.

产品名称: Phospho-JNK1/2/3(T183+T183+T221) Rabbit
Monoclonal Antibody
产品货号: AMRe83788



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunohistochemical analysis of paraffin-embedded human breast, using Phospho-JNK1/2/3 (T183+T183+T221) Antibody.

注意事项 (Note)

For research use only .