

产品名称: EEA1 Rabbit Monoclonal Antibody
产品货号: AMRe83792



产品概述 (Summary)

产品名称 (Production Name)	EEA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEA1
别名 (Alternative Names)	Early endosome antigen 1; EEA1; MST105; ZFYVE2;; Early endosome antigen 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15075.A synthesized peptide derived from human Early endosome antigen 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 162 kDa ; Observed MW: 170 kDa

研究背景 (Background)


Binds phospholipid vesicles containing phosphatidylinositol 3-phosphate and participates in endosomal trafficking.

产品名称: EEA1 Rabbit Monoclonal Antibody
产品货号: AMRe83792




研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: EEA1 Rabbit Monoclonal Antibody
产品货号: AMRe83792



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .