

产品名称: NALP3 Rabbit Monoclonal Antibody
产品货号: AMRe83802



产品概述 (Summary)

产品名称 (Production Name)	NALP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NALP3
别名 (Alternative Names)	FCU; MWS; FCAS; Cias1; Mmig1; NLRP3; Pypaf1; Ail/AVP; AGTAVPRL; Cryopyrin;;NLRP3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96P20.A synthesized peptide derived from human NLRP3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	118 kDa

研究背景 (Background)


产品名称: NALP3 Rabbit Monoclonal Antibody
产品货号: AMRe83802




Sensor component of the NLRP3 inflammasome, which mediates inflammasome activation in response to defects in membrane integrity, leading to secretion of inflammatory cytokines IL1B and IL18 and pyroptosis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: NALP3 Rabbit Monoclonal Antibody
产品货号: AMRe83802




Kynurenic acid/GPR35 axis restricts NLRP3 inflammasome activation and exacerbates colitis in mice with social stress. -BRAIN

BEHAVIOR AND IMMUNITY



Immunofluorescent analysis of Raw264.7 cells, using NALP3 Antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma, using NALP3 Antibody.

注意事项 (Note)

For research use only .