

产品名称: Integrin alpha E Rabbit Monoclonal Antibody
产品货号: AMRe83804



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha E Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin alpha E
别名 (Alternative Names)	Integrin alpha-E; HML-1 antigen; Integrin alpha-IEL; Mucosal lymphocyte 1 antigen; CD103; ITGAE;Integrin alpha E
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P38570.A synthesized peptide derived from human Integrin alpha E

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	130 kDa

研究背景 (Background)


产品名称: Integrin alpha E Rabbit Monoclonal Antibody
产品货号: AMRe83804



Integrin alpha-E/beta-7 is a receptor for E-cadherin. It mediates adhesion of intra-epithelial T-lymphocytes to epithelial cell monolayers.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil, using Integrin alpha E Antibody.

注意事项 (Note)

For research use only .