

产品名称: FNTB Rabbit Monoclonal Antibody
产品货号: AMRe83805



产品概述 (Summary)

产品名称 (Production Name)	FNTB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FNTB CAAX farnesyltransferase beta subunit; EC 2.5.1.58; farnesyltransferase beta;
别名 (Alternative Names)	FNTB; PFTB; Protein farnesyltransferase beta subunit; RAS proteins prenyltransferase beta;; FNTB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49356.A synthesized peptide derived from human FNTB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa ; Observed MW: 46 kDa

产品名称: FNTB Rabbit Monoclonal Antibody
产品货号: AMRe83805




研究背景 (Background)


Essential subunit of the farnesyltransferase complex. Catalyzes the transfer of a farnesyl moiety from farnesyl diphosphate to a cysteine at the fourth position from the C-terminus of several proteins having the C-terminal sequence Cys-aliphatic-aliphatic-X.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of FNTB expression in (1) HepG2 cell lysate; (2) K562 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



产品名称: FNTB Rabbit Monoclonal Antibody
产品货号: AMRe83805



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .