

产品名称: hnRNP C1/C2 Rabbit Monoclonal Antibody
产品货号: AMRe83809



产品概述 (Summary)

产品名称 (Production Name)	hnRNP C1/C2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	hnRNP C1/C2
别名 (Alternative Names)	HNRNPC; HNRPC; Heterogeneous nuclear ribonucleoproteins C1/C2; ROC; hnRNP C1 / hnRNP C2; hnRNP-C1/C2;;hnRNP C1/C2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P07910.A synthesized peptide derived from human hnRNP C1/C2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa ; Observed MW: 40 kDa

产品名称: hnRNP C1/C2 Rabbit Monoclonal Antibody
产品货号: AMRe83809




研究背景 (Background)


Binds pre-mRNA and nucleates the assembly of 40S hnRNP particles. Interacts with poly-U tracts in the 3'-UTR or 5'-UTR of mRNA and modulates the stability and the level of translation of bound mRNA molecules.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of hnRNP C1/C2 expression in (1) HeLa cell lysate; (2) HepG2 cell lysate.



All lanes use the Antibody at 1:5000 dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded human colon, using hnRNP C1/C2 Antibody.

产品名称: hnRNP C1/C2 Rabbit Monoclonal Antibody
产品货号: AMRe83809



注意事项 (Note)

For research use only .