

产品名称: HMBS Rabbit Monoclonal Antibody
产品货号: AMRe83811



产品概述 (Summary)

产品名称 (Production Name)	HMBS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMBS
别名 (Alternative Names)	UPS; PBGD; PORC; PBG-D; HMBS; Hydroxymethylbilane synthase;;Porphobilinogen deaminase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P08397.A synthesized peptide derived from human Porphobilinogen deaminase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa ; Observed MW: 42 kDa

产品名称: HMBS Rabbit Monoclonal Antibody
产品货号: AMRe83811




研究背景 (Background)


As part of the heme biosynthetic pathway, catalyzes the sequential polymerization of four molecules of porphobilinogen to form hydroxymethylbilane, also known as preuroporphyrinogen.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of HMBS expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: HMBS Rabbit Monoclonal Antibody
产品货号: AMRe83811



注意事项 (Note)

For research use only .