

产品名称: MAD2L1 binding protein Rabbit Monoclonal Antibody
产品货号: AMRe83812



产品概述 (Summary)

产品名称 (Production Name)	MAD2L1 binding protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAD2L1 binding protein
别名 (Alternative Names)	Caught by MAD2 protein; CMT2; KIAA0110; MAD2L1 binding protein; MAD2L1-binding protein; MAD2L1BP; MD2BP; MGC11282; RP1-261G23.6;;MAD2L1BP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15013.A synthesized peptide derived from human MAD2L1BP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa ; Observed MW: 34 kDa

产品名称: MAD2L1 binding protein Rabbit Monoclonal Antibody
产品货号: AMRe83812




研究背景 (Background)


May function to silence the spindle checkpoint and allow mitosis to proceed through anaphase by binding MAD2L1 after it has become dissociated from the MAD2L1-CDC20 complex.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MAD2L1BP expression in A431 cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using MAD2L1BP Antibody.

注意事项 (Note)

For research use only .