

产品名称: GART Rabbit Monoclonal Antibody
产品货号: AMRe83813



产品概述 (Summary)

产品名称 (Production Name)	GART Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GART
别名 (Alternative Names)	GARS; GART; AIRS; PRGS; Trifunctional purine biosynthetic protein adenosine 3;;GART
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P22102.A synthesized peptide derived from human GART

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 108 kDa ; Observed MW: 107 kDa

产品名称: GART Rabbit Monoclonal Antibody
产品货号: AMRe83813




研究背景 (Background)


Trifunctional enzyme that catalyzes three distinct reactions as part of the 'de novo' inosine monophosphate biosynthetic pathway.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of GART expression in HeLa cell lysate.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: GART Rabbit Monoclonal Antibody
产品货号: AMRe83813



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded human cervix cancer, using GART Antibody.

注意事项 (Note)

For research use only .