

产品名称: WASF2 Rabbit Monoclonal Antibody
产品货号: AMRe83814



产品概述 (Summary)

产品名称 (Production Name)	WASF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WASF2 Wiskott-Aldrich syndrome protein family member 2WASP family protein member 2; Protein WAVE-2; Verprolin homology domain-containing protein 2; WASF2;;WASF2
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y6W5.A synthesized peptide derived from human WASF2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa ; Observed MW: 80 kDa

产品名称: WASF2 Rabbit Monoclonal Antibody
产品货号: AMRe83814




研究背景 (Background)


Downstream effector molecule involved in the transmission of signals from tyrosine kinase receptors and small GTPases to the actin cytoskeleton. Promotes formation of actin filaments. Part of the WAVE complex that regulates lamellipodia formation. The WAVE complex regulates actin filament reorganization via its interaction with the Arp2/3 complex.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




Western blot analysis of WASF2 expression in K562 cell lysate.



产品名称: WASF2 Rabbit Monoclonal Antibody
产品货号: AMRe83814



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunohistochemical analysis of paraffin-embedded human colon, using WASF2 Antibody.

注意事项 (Note)

For research use only .