

**产品名称: NSD3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe83815**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | NSD3 Rabbit Monoclonal Antibody        |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,FC                                  |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                     |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | NSD3                                                                                                                                |
| 别名 (Alternative Names) | Wolf Hirschhorn syndrome candidate 1 like 1; WHSC1L1; NSD3; pp14328; Histone-lysine N-methyltransferase NSD3; Protein whistle;;NSD3 |
| 基因 ID (Gene ID)        |                                                                                                                                     |
| 蛋白 ID (SwissProt ID)   | Q9BZ95.A synthesized peptide derived from human NSD3                                                                                |

## 产品应用 (Application)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000,FC 1:20-1:100                |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 162 kDa ; Observed MW: 161 kDa |

## 研究背景 (Background)


**产品名称: NSD3 Rabbit Monoclonal Antibody**  
**产品货号: AMRe83815**



Histone methyltransferase. Preferentially dimethylates 'Lys-4' and 'Lys-27' of histone H3 forming H3K2me2 and H3K27me2. H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation, while 'Lys-27' is a mark for transcriptional repression.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of NSD3 expression in MCF-7 cell lysate.

## 注意事项 (Note)

For research use only .