

产品名称: KCNA1 Rabbit Monoclonal Antibody
产品货号: AMRe83816



产品概述 (Summary)

产品名称 (Production Name)	KCNA1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KCNA1
别名 (Alternative Names)	HBK1; HUK1; EA1; MK1; AEMK; MBK1; RBK1; KCNA1; potassium channel Kv1.1; potassium voltage-gated channel subfamily A member 1;;KCNA1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q09470.A synthesized peptide derived from human KCNA1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	56 kDa

研究背景 (Background)


产品名称: KCNA1 Rabbit Monoclonal Antibody
产品货号: AMRe83816



Voltage-gated potassium channel that mediates transmembrane potassium transport in excitable membranes, primarily in the brain and the central nervous system, but also in the kidney .Contributes to the regulation of the membrane potential and nerve signaling, and prevents neuronal hyperexcitability.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of KCNA1 expression in Human fetal brain lysate.

注意事项 (Note)

For research use only .