

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe83821



产品概述 (Summary)

产品名称 (Production Name)	Hsp70 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp70 Heat shock 70 kDa protein 1; Heat shock 70 kDa protein 1/2; Heat shock 70 kDa protein 1A; Heat shock 70 kDa protein 1A/1B; Heat shock 70kDa protein 1B; Heat shock induced protein; HEL S 103; HSP70 1; HSP70 1A; HSP70 1B; HSP70 2; HSP70-1; HSP70-1/HSP70-2; HSP70-1A; HSP70-2; HSP70.1; HSP70.1/HSP70.2; HSP72; HSPA1; HSPA1A; HSPA1B;;Hsp70
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P0DMV8/P0DMV9.A synthesized peptide derived from human Hsp70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-
----------------------	--

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe83821



1:100, IP 1:20-1:50


蛋白分子量 (Molecular Weight) 70 kDa

研究背景 (Background)


Molecular chaperone implicated in a wide variety of cellular processes, including protection of the proteome from stress, folding and transport of newly synthesized polypeptides, activation of proteolysis of misfolded proteins and the formation and dissociation of protein complexes.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Hsp70 expression in (1) HeLa cell lysate; (2) MCF7 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: Hsp70 Rabbit Monoclonal Antibody
产品货号: AMRe83821



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)

注意事项 (Note)

For research use only .