

产品名称: CD36 Rabbit Monoclonal Antibody
产品货号: AMRe83833



产品概述 (Summary)

产品名称 (Production Name)	CD36 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD36 Adipocyte membrane protein; CD36; CD36; CD36 antigen (collagen type I receptor, thrombospondin receptor); CD36 antigen; CD36 molecule (thrombospondin receptor); CD36 molecule; CD36_HUMAN; CHDS7; Cluster determinant 36; Collagen receptor;;CD36
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P16671.A synthesized peptide derived from human CD36

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa ; Observed MW: 70-120 kDa

产品名称: CD36 Rabbit Monoclonal Antibody
产品货号: AMRe83833




研究背景 (Background)


Multifunctional glycoprotein that acts as a receptor for a broad range of ligands. Ligands can be of proteinaceous nature like thrombospondin, fibronectin, collagen or amyloid-beta as well as of lipidic nature such as oxidized low-density lipoprotein (oxLDL), anionic phospholipids, long-chain fatty acids and bacterial diacylated lipopeptides.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of CD36 expression in 3T3 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: CD36 Rabbit Monoclonal Antibody
产品货号: AMRe83833



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .