

产品名称: IRF6 Rabbit Monoclonal Antibody
产品货号: AMRe83836



产品概述 (Summary)

产品名称 (Production Name)	IRF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF6
别名 (Alternative Names)	Interferon regulatory factor 6; IRF 6; LPS;O6; PIT; PPS; PPS1; VWS; VWS1;;IRF6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14896.A synthesized peptide derived from human IRF6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	53 kDa

研究背景 (Background)

Probable DNA-binding transcriptional activator. Key determinant of the keratinocyte proliferation-differentiation switch


产品名称: IRF6 Rabbit Monoclonal Antibody
产品货号: AMRe83836




involved in appropriate epidermal development.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of IRF6 expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: IRF6 Rabbit Monoclonal Antibody
产品货号: AMRe83836



注意事项 (Note)

For research use only .