

产品名称: Cdk8 Rabbit Monoclonal Antibody
产品货号: AMRe83837



产品概述 (Summary)

产品名称 (Production Name)	Cdk8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cdk8 CDK8 protein kinase; Cell division protein kinase 8; Cyclin Dependent kinase 8;
别名 (Alternative Names)	K35; Mediator complex subunit cdk8; Mediator of RNA polymerase II transcription subunit cdk8; Protein kinase K35;;CDK8
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49336.A synthesized peptide derived from human CDK8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	53 kDa

产品名称: Cdk8 Rabbit Monoclonal Antibody
产品货号: AMRe83837




研究背景 (Background)


Component of the Mediator complex, a coactivator involved in regulated gene transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




产品名称: Cdk8 Rabbit Monoclonal Antibody
产品货号: AMRe83837




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:500 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .