

产品名称: 14-3-3 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83841



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 gamma Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	14-3-3 gamma
别名 (Alternative Names)	1433G, 143G, KCIP-1, Protein kinase C inhibitor protein-1, YWHAG; gamma polypeptide;;14-3-3 gamma
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61981.A synthesized peptide derived from human 14-3-3 gamma

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	28 kDa

研究背景 (Background)


产品名称: 14-3-3 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83841




Induce target protein conformational changes that modify target protein function. Distinct temporal and spatial expression patterns of 14-3-3 isoforms have been observed in development and in acute response to extracellular signals and drugs, suggesting that 14-3-3 isoforms may perform different functions despite their sequence similarities.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of 14-3-3 gamma expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: 14-3-3 gamma Rabbit Monoclonal Antibody
产品货号: AMRe83841



注意事项 (Note)

For research use only .