

产品名称: ATF7 Rabbit Monoclonal Antibody
产品货号: AMRe83842



产品概述 (Summary)

产品名称 (Production Name)	ATF7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATF7 Cyclic AMP-dependent transcription factor ATF-7; cAMP-dependent
别名 (Alternative Names)	transcription factor ATF-7; Activating transcription factor 7; Transcription factor ATF-A; ATF7; ATFA;;ATF7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P17544.A synthesized peptide derived from human ATF7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa ; Observed MW: 53 kDa

产品名称: ATF7 Rabbit Monoclonal Antibody
产品货号: AMRe83842




研究背景 (Background)


Stress-responsive chromatin regulator that plays a role in various biological processes including innate immunological memory, adipocyte differentiation or telomerase regulation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of ATF7 expression in (1) Raji cell lysate; (2) Raw 264.7 cell lysate.

产品名称: ATF7 Rabbit Monoclonal Antibody
产品货号: AMRe83842



注意事项 (Note)

For research use only .