

产品名称: Phospho-JunD(S255) Rabbit Monoclonal Antibody
产品货号: AMRe83848



产品概述 (Summary)

产品名称 (Production Name)	Phospho-JunD(S255) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-JunD(S255)
别名 (Alternative Names)	Activator protein 1; AP 1; Jun D; jun D proto oncogene; Jund; JunD FL isoform; Transcription factor jun D;;p-JunD (S255)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P17535.A synthesized peptide derived from human JunD around the phosphorylation site of S255

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa ; Observed MW: 38,42 kDa

产品名称: Phospho-JunD(S255) Rabbit Monoclonal Antibody
产品货号: AMRe83848




研究背景 (Background)


Transcription factor binding AP-1 sites. Heterodimerizes with proteins of the FOS family to form an AP-1 transcription factor complex, thereby enhancing their DNA binding activity to an AP-1 consensus sequence 3'-TGA[GC]TCA-5' and enhancing their transcriptional activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-JunD (S255) expression in (1) HeLa cell lysate treated with Alkaline Phosphatase; (2) HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .