

产品名称: HDAC8 Rabbit Monoclonal Antibody
产品货号: AMRe83858



产品概述 (Summary)

产品名称 (Production Name)	HDAC8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HDAC8
别名 (Alternative Names)	Histone deacetylase 8; HD8; HDAC8; HDACL1; CDA07; RPD 3; MRXS6;;HDAC8
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BY41.A synthesized peptide derived from human HDAC8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	42 kDa

研究背景 (Background)


产品名称: HDAC8 Rabbit Monoclonal Antibody
产品货号: AMRe83858




Histone deacetylase that catalyzes the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of HDAC8 expression in HeLa cell lysate.

产品名称: HDAC8 Rabbit Monoclonal Antibody
产品货号: AMRe83858



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .