

产品名称: HSPA14 Rabbit Monoclonal Antibody
产品货号: AMRe83859



产品概述 (Summary)

产品名称 (Production Name)	HSPA14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA14
别名 (Alternative Names)	HSPA14; Heat shock protein HSP60; Heat shock 70 kDa protein 14; HSP70L1; HSP60; Heat shock 70kDa protein 14; HSP70-like protein 1; HSP70-4;;Hsp70 14
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q0VDF9.A synthesized peptide derived from human Hsp70 14

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	55 kDa

产品名称: HSPA14 Rabbit Monoclonal Antibody
产品货号: AMRe83859




研究背景 (Background)

Component of the ribosome-associated complex (RAC), a complex involved in folding or maintaining nascent polypeptides in a folding-competent state. In the RAC complex, binds to the nascent polypeptide chain, while DNAJC2 stimulates its ATPase activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSPA14 expression in K562 cell lysate.

注意事项 (Note)

For research use only .