

产品名称: SUMO 4 Rabbit Monoclonal Antibody
产品货号: AMRe83861



产品概述 (Summary)

产品名称 (Production Name)	SUMO 4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SUMO 4
别名 (Alternative Names)	IDDM5; SMT3H4; SUMO-4; dJ281H8.4;; SUMO 4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q6EEV6.A synthesized peptide derived from human SUMO 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 11 kDa ; Observed MW: 15 kDa

研究背景 (Background)

Ubiquitin-like protein which can be covalently attached to target lysines as a monomer. Does not seem to be involved in


产品名称: SUMO 4 Rabbit Monoclonal Antibody
产品货号: AMRe83861




protein degradation and may modulate protein subcellular localization, stability or activity. Upon oxidative stress, conjugates to various anti-oxidant enzymes, chaperones, and stress defense proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SUMO4 expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .