

产品名称: CDC5 Rabbit Monoclonal Antibody
产品货号: AMRe83863



产品概述 (Summary)

产品名称 (Production Name)	CDC5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDC5
别名 (Alternative Names)	CDC5; CEF1; PCDC5RP; CDC5-LIKE; dJ319D22.1; CDC5L;;CDC5L
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99459.A synthesized peptide derived from human CDC5L

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa ; Observed MW: 100 kDa

研究背景 (Background)

DNA-binding protein involved in cell cycle control. May act as a transcription activator. Plays a role in pre-mRNA splicing as


产品名称: CDC5 Rabbit Monoclonal Antibody
产品货号: AMRe83863




core component of precatalytic, catalytic and postcatalytic spliceosomal complexes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CDC5 expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .