

产品名称: Cdc27/APC3 Rabbit Monoclonal Antibody
产品货号: AMRe83864



产品概述 (Summary)

产品名称 (Production Name)	Cdc27/APC3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cdc27/APC3
别名 (Alternative Names)	CDC27Hs; Cell division cycle protein 27 homolog; HNUC; APC3; ANAPC3;;CDC27
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P30260.A synthesized peptide derived from human CDC27

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa ; Observed MW: 91 kDa

研究背景 (Background)


产品名称: Cdc27/APC3 Rabbit Monoclonal Antibody
产品货号: AMRe83864




Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. The APC/C complex acts by mediating ubiquitination and subsequent degradation of target proteins: it mainly mediates the formation of 'Lys-11'-linked polyubiquitin chains and, to a lower extent, the formation of 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cdc27/APC3 expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .