

产品名称: Cdc23/APC8 Rabbit Monoclonal Antibody
产品货号: AMRe83868



产品概述 (Summary)

产品名称 (Production Name)	Cdc23/APC8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cdc23/APC8
别名 (Alternative Names)	CDC23; ANAPC8; APC8; Cell division cycle 23; CUT23; Cyclosome subunit 8;; CDC23
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UJX2.A synthesized peptide derived from human CDC23

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa ; Observed MW: 64 kDa

研究背景 (Background)


产品名称: Cdc23/APC8 Rabbit Monoclonal Antibody
产品货号: AMRe83868




Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. The APC/C complex acts by mediating ubiquitination and subsequent degradation of target proteins: it mainly mediates the formation of 'Lys-11'-linked polyubiquitin chains and, to a lower extent, the formation of 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Cdc23/APC8 expression in (1)HepG2 cell lysate; (2)Jurkat cell lysate.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




产品名称: Cdc23/APC8 Rabbit Monoclonal Antibody
产品货号: AMRe83868



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .