

产品名称: Hsp20 Rabbit Monoclonal Antibody
产品货号: AMRe83870



产品概述 (Summary)

产品名称 (Production Name)	Hsp20 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp20
别名 (Alternative Names)	HSPB6; HEL55; Hsp20; Heat-shock 20 kDa like-protein p20;;Hsp20
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14558.A synthesized peptide derived from human Hsp20

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)

Small heat shock protein which functions as a molecular chaperone probably maintaining denatured proteins in a folding-


产品名称: Hsp20 Rabbit Monoclonal Antibody
产品货号: AMRe83870



competent state. Seems to have versatile functions in various biological processes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Hsp20 expression in mouse skeletal muscle lysate.

注意事项 (Note)

For research use only .