

产品名称: Hsp47 Rabbit Monoclonal Antibody
产品货号: AMRe83873



产品概述 (Summary)

产品名称 (Production Name)	Hsp47 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp47
别名 (Alternative Names)	CBP1; CBP2; OI10; gp46; AsTP3; HSP47; PIG14; PPRM; RA-A47; SERPINH2;;Serp H1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50454.A synthesized peptide derived from human Serpin H1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	46 kDa

研究背景 (Background)


产品名称: Hsp47 Rabbit Monoclonal Antibody
产品货号: AMRe83873




This gene encodes a member of the serpin superfamily of serine proteinase inhibitors. The encoded protein is localized to the endoplasmic reticulum and plays a role in collagen biosynthesis as a collagen-specific molecular chaperone. Autoantibodies to the encoded protein have been found in patients with rheumatoid arthritis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: Hsp47 Rabbit Monoclonal Antibody
产品货号: AMRe83873



注意事项 (Note)

For research use only .