

产品名称: delta 2 Catenin Rabbit Monoclonal Antibody
产品货号: AMRe83879



产品概述 (Summary)

产品名称 (Production Name)	delta 2 Catenin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	delta 2 Catenin ARM-repeat protein; Catenin delta 2; CTNND 2; Ctnnd2; Delta catenin; Delta-catenin; GT24; Neural plakophilin related arm repeat protein; Neurojuning; NPRAP; T cell delta catenin;;delta 2 Catenin
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O35927(mouse)(human:Q9UQB3).A synthesized peptide derived from mouse delta 2 Catenin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	135 kDa

产品名称: delta 2 Catenin Rabbit Monoclonal Antibody
产品货号: AMRe83879




研究背景 (Background)


Has a critical role in neuronal development, particularly in the formation and/or maintenance of dendritic spines and synapses. Involved in the regulation of canonical Wnt signaling. It probably acts on beta-catenin turnover, facilitating beta-catenin interaction with GSK3B, phosphorylation, ubiquitination and degradation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of delta 2 Catenin expression in Mouse brain cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .