

产品名称: EBP50 Rabbit Monoclonal Antibody
产品货号: AMRe83885



产品概述 (Summary)

产品名称 (Production Name)	EBP50 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EBP50
别名 (Alternative Names)	SLC9A3R1; EBP50; NPHLOP2; NHERF; NHERF-1; NHERF1;;NHERF 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14745.A synthesized peptide derived from human NHERF 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa ; Observed MW: 50 kDa

研究背景 (Background)

Scaffold protein that connects plasma membrane proteins with members of the ezrin/moesin/radixin family and thereby


产品名称: EBP50 Rabbit Monoclonal Antibody
产品货号: AMRe83885




helps to link them to the actin cytoskeleton and to regulate their surface expression. Necessary for recycling of internalized ADRB2.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of EBP50 expression in (1) HepG2 cell lysate; (2) PC12 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: EBP50 Rabbit Monoclonal Antibody
产品货号: AMRe83885



注意事项 (Note)

For research use only .