

产品名称: PGC1alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe83890



产品概述 (Summary)

产品名称 (Production Name)	PGC1alpha/beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGC1alpha/beta
别名 (Alternative Names)	PPARGC1A; LEM6; PGC1; PGC1A; PGC-1v; PPARGC1; PGC-1(alpha); PGC1 alpha/beta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBK2.A synthesized peptide derived from human PGC1 alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 91,113 kDa ; Observed MW: 113 kDa

研究背景 (Background)


产品名称: PGC1alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe83890




Transcriptional coactivator for steroid receptors and nuclear receptors. Greatly increases the transcriptional activity of PPARG and thyroid hormone receptor on the uncoupling protein promoter. Can regulate key mitochondrial genes that contribute to the program of adaptive thermogenesis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1000 dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1000 dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1000 dilution for 1 hour at room temperature.

产品名称: PGC1alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe83890



Immunofluorescent analysis of HeLa cells using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .