

产品名称: PDIA6 Rabbit Monoclonal Antibody
产品货号: AMRe83891



产品概述 (Summary)

产品名称 (Production Name)	PDIA6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDIA6
别名 (Alternative Names)	PDIA6; ERP5, P5, TXNDC7; ER protein 5; ERp5;; PDIA6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15084.A synthesized peptide derived from human PDIA6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	48 kDa

研究背景 (Background)

May function as a chaperone that inhibits aggregation of misfolded proteins. Negatively regulates the unfolded protein


产品名称: PDIA6 Rabbit Monoclonal Antibody
产品货号: AMRe83891




response (UPR) through binding to UPR sensors such as ERN1, which in turn inactivates ERN1 signaling. May also regulate the UPR via the EIF2AK3 UPR sensor.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of PDIA6 expression in K562 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: PDIA6 Rabbit Monoclonal Antibody
产品货号: AMRe83891



注意事项 (Note)

For research use only .