

产品名称: LYVE1 Rabbit Monoclonal Antibody
产品货号: AMRe83893



产品概述 (Summary)

产品名称 (Production Name)	LYVE1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LYVE1
别名 (Alternative Names)	CRSBP 1; CRSBP-1; CRSBP1; hyaluronic acid receptor;;LYVE1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y5Y7.A synthesized peptide derived from human LYVE1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)

Ligand-specific transporter trafficking between intracellular organelles (TGN) and the plasma membrane. Plays a role in


产品名称: LYVE1 Rabbit Monoclonal Antibody
产品货号: AMRe83893




autocrine regulation of cell growth mediated by growth regulators containing cell surface retention sequence binding (CRS). May act as a hyaluronan (HA) transporter, either mediating its uptake for catabolism within lymphatic endothelial cells themselves, or its transport into the lumen of afferent lymphatic vessels for subsequent re-uptake and degradation in lymph nodes.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: LYVE1 Rabbit Monoclonal Antibody
产品货号: AMRe83893



Western blot analysis of LYVE1 expression in A549 cell lysate.

注意事项 (Note)

For research use only .