

产品名称: 58K Golgi protein Rabbit Monoclonal Antibody
产品货号: AMRe83895



产品概述 (Summary)

产品名称 (Production Name)	58K Golgi protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	58K Golgi protein
别名 (Alternative Names)	Formimidoyltetrahydrofolate cyclodeaminase; FTCD; Glutamate formiminotransferase; LCHC 1;;FTCD
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95954.A synthesized peptide derived from human FTCD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa ; Observed MW: 58 kDa

研究背景 (Background)


产品名称: 58K Golgi protein Rabbit Monoclonal Antibody
产品货号: AMRe83895



Folate-dependent enzyme, that displays both transferase and deaminase activity. Serves to channel one-carbon units from formiminoglutamate to the folate pool. Binds and promotes bundling of vimentin filaments originating from the Golgi.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 58K Golgi protein expression in human fetal liver lysate.

注意事项 (Note)

For research use only .