

产品名称: Dystrophin Rabbit Monoclonal Antibody
产品货号: AMRe83900



产品概述 (Summary)

产品名称 (Production Name)	Dystrophin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Dystrophin
别名 (Alternative Names)	BMD; CMD3B; DMD; Dystrophin;;Dystrophin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P11532.A synthesized peptide derived from human Dystrophin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	427 kDa

研究背景 (Background)

Anchors the extracellular matrix to the cytoskeleton via F-actin. Ligand for dystroglycan. Component of the dystrophin-


产品名称: Dystrophin Rabbit Monoclonal Antibody
产品货号: AMRe83900



associated glycoprotein complex which accumulates at the neuromuscular junction (NMJ) and at a variety of synapses in the peripheral and central nervous systems and has a structural function in stabilizing the sarcolemma. Also implicated in signaling events and synaptic transmission.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Dystrophin expression in human fetal heart lysate.

注意事项 (Note)

For research use only .