

产品名称: HLA G Rabbit Monoclonal Antibody
产品货号: AMRe83901



产品概述 (Summary)

产品名称 (Production Name)	HLA G Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HLA G
别名 (Alternative Names)	HLA G antigen; MHC class I antigen G; HLA-G; HLA60; HLAG;;HLA G
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P17693.A synthesized peptide derived from human HLA G

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa ; Observed MW: 42 kDa

研究背景 (Background)

Non-classical major histocompatibility class Ib molecule involved in immune regulatory processes at the maternal-fetal


产品名称: HLA G Rabbit Monoclonal Antibody
产品货号: AMRe83901




interface. In complex with B2M/beta-2 microglobulin binds a limited repertoire of nonamer self-peptides derived from intracellular proteins including histones and ribosomal proteins.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of HLA G expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: HLA G Rabbit Monoclonal Antibody
产品货号: AMRe83901



注意事项 (Note)

For research use only .