

产品名称: NUMB Rabbit Monoclonal Antibody
产品货号: AMRe83902



产品概述 (Summary)

产品名称 (Production Name)	NUMB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUMB
别名 (Alternative Names)	c14_5527; h Numb; Numb; Numb homolog; S171; ;NUMB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P49757.A synthesized peptide derived from human NUMB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	71 kDa

研究背景 (Background)


产品名称: NUMB Rabbit Monoclonal Antibody
产品货号: AMRe83902




Numb acts as a negative regulator of Notch signaling by promoting ubiquitination and degradation of Notch. The protein is asymmetrically segregated into one daughter cell during cell division, producing two daughter cells with different responses to Notch signaling and different cell fates. The localization of Numb can also be regulated by G-protein coupled receptor (GPCR) and PKC signaling.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of NUMB expression in Ramos cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



产品名称: NUMB Rabbit Monoclonal Antibody
产品货号: AMRe83902



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .