

产品名称: TGN46 Rabbit Monoclonal Antibody
产品货号: AMRe83903



产品概述 (Summary)

产品名称 (Production Name)	TGN46 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TGN46
别名 (Alternative Names)	TGN38; TGN46; TGN48; TGN51; TGOLN2; TTN2; TGN46
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43493.A synthesized peptide derived from human TGN46

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa ; Observed MW: 90 kDa

研究背景 (Background)


May be involved in regulating membrane traffic to and from trans-Golgi network.

产品名称: TGN46 Rabbit Monoclonal Antibody
产品货号: AMRe83903




研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of TGN46 expression in HeLa cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

注意事项 (Note)

产品名称: TGN46 Rabbit Monoclonal Antibody
产品货号: AMRe83903



For research use only .