

产品名称: CNPase Rabbit Monoclonal Antibody
产品货号: AMRe83904



产品概述 (Summary)

产品名称 (Production Name)	CNPase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CNPase
别名 (Alternative Names)	CNP; CNP1; CNPase;;CNPase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P09543.A synthesized peptide derived from human CNPase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	48 kDa

研究背景 (Background)

Catalyzes the formation of 2'-nucleotide products from 2',3'-cyclic substrates. May participate in RNA metabolism in the


产品名称: CNPase Rabbit Monoclonal Antibody
产品货号: AMRe83904




myelinating cell, CNP is the third most abundant protein in central nervous system myelin.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CNPase expression in C6 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .