

产品名称: Endothelin B Receptor Rabbit Monoclonal Antibody
产品货号: AMRe83912



产品概述 (Summary)

产品名称 (Production Name)	Endothelin B Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Endothelin B Receptor
别名 (Alternative Names)	ABCD5; Ednra; EDNRB; Endothelin B receptor precursor; ETB; ETBR; ETRB; Hirschsprung disease;;EDNRB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P24530.A synthesized peptide derived from human EDNRB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	50 kDa

产品名称: Endothelin B Receptor Rabbit Monoclonal Antibody
产品货号: AMRe83912




研究背景 (Background)

Non-specific receptor for endothelin 1, 2, and 3. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Endothelin B Receptor expression in 293T cell lysate.

注意事项 (Note)

For research use only .