

产品名称: FABP1 Rabbit Monoclonal Antibody
产品货号: AMRe83914



产品概述 (Summary)

产品名称 (Production Name)	FABP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FABP1
别名 (Alternative Names)	FABP1; FABPL; Fatty acid binding protein 1 liver; L FABP;;FABP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P07148.A synthesized peptide derived from human FABP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	14 kDa

研究背景 (Background)

Plays a role in lipoprotein-mediated cholesterol uptake in hepatocytes. Binds cholesterol. Binds free fatty acids and their


产品名称: FABP1 Rabbit Monoclonal Antibody
产品货号: AMRe83914




coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of FABP1 expression in human fetal liver lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: FABP1 Rabbit Monoclonal Antibody
产品货号: AMRe83914



注意事项 (Note)

For research use only .