

产品名称: CYP2E1 Rabbit Monoclonal Antibody
产品货号: AMRe83918



产品概述 (Summary)

产品名称 (Production Name)	CYP2E1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYP2E1
别名 (Alternative Names)	Cytochrome P450 2E1; CYP2E1; Cytochrome P450-J; CYP2E; P450C2E; CYP2E1;;Cytochrome P450 2E1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05181.A synthesized peptide derived from human Cytochrome P450 2E1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	57 kDa

研究背景 (Background)


产品名称: CYP2E1 Rabbit Monoclonal Antibody
产品货号: AMRe83918




A cytochrome P450 monooxygenase involved in the metabolism of fatty acids. Mechanistically, uses molecular oxygen inserting one oxygen atom into a substrate, and reducing the second into a water molecule, with two electrons provided by NADPH via cytochrome P450 reductase (NADPH--hemoprotein reductase).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CYP2E1 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .