

产品名称: TIMP2 Rabbit Monoclonal Antibody
产品货号: AMRe83919



产品概述 (Summary)

产品名称 (Production Name)	TIMP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIMP2
别名 (Alternative Names)	CSC21K; TIMP 2; TIMP2;;TIMP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P16035.A synthesized peptide derived from human TIMP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa ; Observed MW: 22 kDa

研究背景 (Background)

Complexes with metalloproteinases (such as collagenases) and irreversibly inactivates them by binding to their catalytic


产品名称: TIMP2 Rabbit Monoclonal Antibody
产品货号: AMRe83919



zinc cofactor. Known to act on MMP-1, MMP-2, MMP-3, MMP-7, MMP-8, MMP-9, MMP-10, MMP-13, MMP-14, MMP-15, MMP-16 and MMP-19.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TIMP2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .