

产品名称: TLR5 Rabbit Monoclonal Antibody
产品货号: AMRe83920



产品概述 (Summary)

产品名称 (Production Name)	TLR5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TLR5
别名 (Alternative Names)	SLEB1; TIL3; Tlr5; Toll like receptor 5 precursor;; TLR5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O60602.A synthesized peptide derived from human TLR5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 98 kDa ; Observed MW: 97 kDa

研究背景 (Background)

Pattern recognition receptor (PRR) located on the cell surface that participates in the activation of innate immunity and


产品名称: TLR5 Rabbit Monoclonal Antibody
产品货号: AMRe83920




inflammatory response. Recognizes small molecular motifs named pathogen-associated molecular pattern (PAMPs) expressed by pathogens and microbe-associated molecular patterns (MAMPs) usually expressed by resident microbiota.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: TLR5 Rabbit Monoclonal Antibody
产品货号: AMRe83920



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:200 dilution.

注意事项 (Note)

For research use only .