

产品名称: PRDX6 Rabbit Monoclonal Antibody
产品货号: AMRe83924



产品概述 (Summary)

产品名称 (Production Name)	PRDX6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRDX6
别名 (Alternative Names)	Peroxiredoxin-6; Antioxidant protein 2; Liver 2D page spot 40; NSGPx; PRDX6; AOP2; KIAA;;PRDX6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P30041.A synthesized peptide derived from human PRDX6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	25 kDa

研究背景 (Background)


产品名称: PRDX6 Rabbit Monoclonal Antibody
产品货号: AMRe83924




Involved in redox regulation of the cell. Can reduce H₂O₂ and short chain organic, fatty acid, and phospholipid hydroperoxides. May play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of PRDX6 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: PRDX6 Rabbit Monoclonal Antibody
产品货号: AMRe83924



注意事项 (Note)

For research use only .