

产品名称: LY6A Rabbit Monoclonal Antibody
产品货号: AMRe83930



产品概述 (Summary)

产品名称 (Production Name)	LY6A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LY6A Ly-6A.2/Ly-6E.1; Ly6a; Ly6al; Lymphocyte antigen 6A-2/6E-1; Lymphocyte
别名 (Alternative Names)	antigen Ly 6A.2/Ly 6E.1; SCA-1; Stem cell antigen 1; T cell activating protein; TAP; ;Ly6a
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05533.A synthesized peptide derived from mouse Ly6a

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa ; Observed MW: 18 kDa

产品名称: LY6A Rabbit Monoclonal Antibody
产品货号: AMRe83930




研究背景 (Background)


T-cell activation.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of LY6A expression in mouse kidney lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: LY6A Rabbit Monoclonal Antibody
产品货号: AMRe83930



注意事项 (Note)

For research use only .