

产品名称: Drebrin Rabbit Monoclonal Antibody
产品货号: AMRe83941



产品概述 (Summary)

产品名称 (Production Name)	Drebrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Drebrin
别名 (Alternative Names)	DOS117E; DBN1; Drebrin 1; Drebrin; Drebrin E2; Drebrin1; DrebrinE; ;Drebrin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q16643.A synthesized peptide derived from human Drebrin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 71 kDa ; Observed MW: 120 kDa

研究背景 (Background)

Actin cytoskeleton-organizing protein that plays a role in the formation of cell projections. Required for actin


产品名称: Drebrin Rabbit Monoclonal Antibody
产品货号: AMRe83941




polymerization at immunological synapses (IS) and for the recruitment of the chemokine receptor CXCR4 to IS.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: Drebrin Rabbit Monoclonal Antibody
产品货号: AMRe83941



注意事项 (Note)

For research use only .