

产品名称: Aggrecan Rabbit Monoclonal Antibody
产品货号: AMRe83950



产品概述 (Summary)

产品名称 (Production Name)	Aggrecan Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Aggrecan
别名 (Alternative Names)	ACAN; AGC1; AGCAN; Aggrecan 1; CSPG1; CSPGCP; JSCATE; Large aggregating proteoglycan; mcspg; mgsk16; MSK 16; SEDK;;Aggrecan
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P16112.A synthesized peptide derived from human Aggrecan

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 261 kDa ; Observed MW: 250 kDa

研究背景 (Background)


产品名称: Aggrecan Rabbit Monoclonal Antibody
产品货号: AMRe83950




This proteoglycan is a major component of extracellular matrix of cartilagenous tissues. A major function of this protein is to resist compression in cartilage. It binds avidly to hyaluronic acid via an N-terminal globular region.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded human cartilage tissue, using Aggrecan Antibody.

注意事项 (Note)

For research use only .